
BIOTESLA LASER

Lasertherapy



electronicamilka@emilka.com
+34 91 463 51 08 - fax: +34 91 470 00 12
ELECTRONICA MILKA S.L.
C/ HERMINIO PUERTAS, 66 - 28011 MADRID - SPAIN

MAIN TREATMENTS

Tendinitis, varicose ulcers, arthrosis, peri-arthritis, tendovaginitis, acute posterior ganglionitis, etc...

Acupuncture Treatments

There is disposable a special lens (optional) to can do treatments in acupuncture points. This lens concentrates the laser beam in only one point.

The treatments have the same effectiveness that traditional needles.

This is especially indicated for patients with fear to the needles.

CHARACTERISTICS:

It has an ultralight support with refrigerator with an intensity of 30 mW (with a pick ok 30 W) with support base. The therapy time can be switched between minutes or seconds. This allow to use acupuncture points.

The laser emission can be selected between continuous or pulsated, being its frequencies from 0 to 5000 Hz.

The timer goes from 1 to 99 minutes or seconds.

The generator also has a sensor to check the laser emission.



The **BIOTESLA** equipment is a infrared diodic laser of half/low intensity and with a wave length of 904 nanometers. It pass throught cutaneous and subcutaneous barrier and operate deeply without burns.

It has an ultralight support with refrigerator with an intensity of 30 mW (with a pick ok 30 W) with support base. The therapy time can be switched between minutes or seconds.

This allow to use acupunture points.

The generator also has a sensor to check the laser emission.

TECHNICAL DATA: BIOTESLA LASER

Supply : 230 V / 50 Hz.
Size : 313x224x95 mm.
Support size: 150x22 mm.
Weight. : 3 Kgs.
Support base.



Guarantee: 1 year

Manufacture and own technical service

Tel.: 91 463 51 08 e-mail.: electronicamilka@emilka.com