

---

# BIOTESLA LASER

Veterinary

---



---

**electronicamilka@emilka.com**  
**+34 91 463 51 08 - fax: +34 91 470 00 12**  
**ELECTRONICA MILKA S.L.**  
**C/ HERMINIO PUERTAS, 66 - 28011 MADRID - SPAIN**

## **MAIN TREATMENTS**

Tendinitis, varicose ulcers, arthrosis, peri-arthritis, tendovaginitis, etc...

### **CHARACTERISTICS:**

It has an ultralight support with refrigerator with an intensity of 30 mW (with a pick ok 30 W) with support base. The therapy time can be switched between minutes or seconds. This allow to use acupunture points.

The laser emission can be selected between continuous or pulsated, being its frequencies from 0 to 5000 Hz.

The timer goes from 1 to 99 minutes or seconds.

The generator also has a sensor to check the laser emission.



The **BIOTESLA** equipment is a infrared diodic laser of half/low intensity and with a wave length of 904 nanometers. It pass throught cutaneous and subcutaneous barrier and operate deeply without burns.

It has an ultralight support with refrigerator with an intensity of 30 mW (with a pick ok 30 W) with support base. The therapy time can be switched between minutes or seconds.

This allow to use acupunture points.

**The generator also has a sensor to check the laser emission.**

## **TECHNICAL DATA: BIOTESLA LASER**

Supply : 230 V / 50 Hz.  
Size : 313x224x95 mm.  
Support size: 150x22 mm.  
Weight. : 3 Kgs.  
Support base.



**Guarantee: 1 year**

**Manufacture and own technical service**

Tel.: 91 463 51 08 e-mail.: [electronicamilka@emilka.com](mailto:electronicamilka@emilka.com)